# Inflatable Chair Pad

## Item No: 748140 / K411

### Instructions for Use

v001: 31/10/25

Please read and retain these instructions for future reference

# 3PAGEN

## ENGLISH

### Dimensions

Seat size – 52 x 35cm

Back size – 49 x 41cm

INFLATABLE CHAIR PAD WITH HAND PUMP SO IT CAN BE PACKED AWAY AND USED WHEREVER YOU NEED IT  
RELIEVE YOUR BACK AND LOWER BACK.

IN ORDER NOT TO SUFFER FROM BACK PAIN IT IS ESSENTIAL TO SIT WELL.

THE FORM FITTING CUSHION SUPPORTS YOUR LOWER BACK AND KEEPS YOUR BACK COMFORTABLE TO RELIEVE ALL TENSION.

VERY EASY TO INFLATE, IT SIMPLY ADAPTS TO ANY CHAIR OR ARMCHAIR AND TRANSFORMS INTO A REAL COCOON OF SOFTNESS FOR YOUR SPINE!

Max load 100kg

### INSTRUCTIONS FOR USE:

Step by Step Instructions & Photos

Unpack the components carefully. Please check all parts have been removed from the packaging.

Inflate:

1. Pull the valve until it comes out.

2. Tilt the valve on its side and insert the nozzle in the valve as shown. Make sure you keep the pump at an angle to the valve as shown. The nozzle will fit inside the valve just move it from side to side. Hold the nozzle end and push gently and start to pump.

3. Pump to the desired air level and firmness then tuck the valve in.

**Check regularly if the pad is at your required air level and firmness.**

The pad has 4 holes as marked and can be tied to a chair using string (not provided).

### Deflate:

Pull out the valve, pinch the mouth of valve and squeeze the cushion seat to deflate.

### CARE AND CLEANING INSTRUCTIONS:

Clean with a dry cloth or duster.

## DISPOSAL

- We use recyclable or recycled packaging where possible.
- Please dispose of all packaging, paper, cartons in accordance with your local authority recycling regulations.
- At the end of the product's lifespan please check with your local council authorised household waste recycling centre for disposal.
- Plastics, polybags – this is made from the following recyclable plastic.

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The Enterprise Department Ltd, Trading as Coopers of Stortford.