



68826 3-Kammer Lungentrainer
3-chamber lung trainer
Entraîneur pulmonaire à 3 chambres
tríkomorový plicní trenažér
3-komorový pľúcny trenažér

DE:

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EN:

Thank you for purchasing our 3-chamber lung trainer.
We hope you are satisfied with your purchase.

How to use your lung trainer

1. hold the device in an upright position
2. inhale normally and then close your lips tightly around the mouthpiece at the end of the green tube. (Fig. 1)
3. LOW FLOW RATE - Breathe in so quickly that only the ball in the first chamber is lifted. The ball in the second chamber must remain in place. This position should be held for three seconds or as long as possible, whichever comes first. (Fig. 2)
4. HIGH FLOW RATE - Inhale fast enough to lift the balls in the first and second chambers. Make sure that the third chamber ball remains in the resting position for the entire duration of this exercise. (Fig. 3)
5. BREATHE OUT - Remove the mouthpiece and breathe out normally (Fig. 4) RELAX (repetition) - After each prolonged deep breath, take a moment to rest and breathe normally.

NOTE: By tilting the device forwards, you can make lung training easier for users who find it difficult to lift the ball or balls while holding the device in an upright position.

Cleaning: If you want to clean the hose, wash it with soap and warm water.

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The Enterprise Department Ltd, Trading as Coopers of Stortford.

UK address: 11 Bridge Street, Bishops Stortford

Herts, UK, CM23 2JU

Tel: 0330 331 0300

Customer Services: enquiries@coopersonstortford.co.uk

Web Site: <https://www.coopersonstortford.co.uk>

The Enterprise Department Ltd, Trading as -

NI address: D.S.B. LD, 160 Bd. De Fourmies, 59100 Roubaix, France

Tel: +44 (0) 330 331 0300

Customer Services: enquiries@coopersonstortford.co.uk

Web Site: <https://www.coopersonstortford.co.uk>

DE/FR/CZ/SK :

3Pagen

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Sankt-Joeris-Strasse 16-28

53477 Alsdorf

Tel: +49 (0) 2404/980-204

Customer Services: kontakt@3pagen.de

Web Site: www.3pagen.de