(GB) FULL BODY MASSAGE MAT WITH HEAT

ASSEMBLING INSTRUCTION

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F4 – Upper legs massage

F5 – Lower legs massage

Each Massage area is powered separately. To activate a specific massage area, press the buttons for each of the four areas separately.

You can activate the massage area or area combination simultaneously.

Body Message Area:

Neck

Back

Lumbar

Legs (Thigh & Calf)

To activate the heat, press the “HEAT” button. To deactivate it, press the heating button again.

When the heating mode is turned on, the LED indicator lights up.

When disconnecting, press the “ON/OFF” button and then remove the plug from the wall power outlet.

The massage mat is equipped with an easy - to - use controller that divides the body into 4 zones for each massage.

You can activate a single massage zone or a continuous sequence of zones.

The LED lights indicate the selected zones.

OUTLET POWER

To connect the electrical power for the massage mat, plug the end of the power cord into the inner part of the wall adapter.

Make sure that all control buttons are in the STOP position before plugging in the wall adapter.

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OPERATION OF THE CONTROLLER

Power & Heating:

ON/OFF – Turns the massage unit on or off.

HEAT – Activates the heating function. Corresponding indicator light will turn on when active.

Mode & Intensity:

MODE F1–F5 – Cycles through 5 preset massage programs.

F1 – Massage heads (upper body/shoulder area)

F2 – Back massage

F3 – Lumbar massage (lower back)

F4 – Upper legs massage

F5 – Lower legs massage

INTENSITY – Adjusts massage intensity level.

L – Low

M – Medium

H – High

Target Area Control:

NECK – Activates massage for the neck area.

BACK – Activates massage for the back area.

LUMBAR – Activates massage for the lumbar (lower back) region.

LEG – Activates massage for the leg area.

Indicator Lights:

F1 – F5 – Lights up to indicate which massage mode is currently selected.

L / M / H – Shows current massage intensity.

ON – Indicates system is powered on.

HEAT – Indicates the heating function is activated.

How to Use:

Turn On the device using the ON/OFF button.

Select Mode using the MODE NECK/BACK/LUMBAR/LEG button to choose your desired massage area/program.

Adjust Intensity with the INTENSITY button (Low, Medium, High).

Activate Heating by pressing the HEAT button if desired.

Use NECK, BACK, LUMBAR, or LEG buttons to directly activate or focus on a specific area.

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