Product name: Nail Fungus Cleaner **3PAGEN**

Item No: 708590

User Manuals

Thoroughly clean all nails before using this device.

Trim nails and thin them as much as possible.

Open the device and simply place the device on the affected nail.

After 7 minutes, the device turns off automatically and one treatment is completed.

Use 2 times a day, each time you can use 1-2 rounds, that is 7 minutes or 14 minutes, stick to it for 3 months and you will see improvement.

Safety & After-sales Service

- 1 The device should be used in strict compliance with the safety and operating precautions provided in this manual.
- 2 Please review the manual before initial use
- 3 No side effects
- 4 Direct or reflected laser radiation should be prevented from eye contact.
- 5 Pregnant women should not use the device close to the womb
- 6 Keep the device away from children.
- 7 Children under the age of 18 should be treated under the supervision of an adult.

Warm Tips

The device takes 3-4 hours to be fully charged. When you turn on the product and find that it does not light up, please charge it first; When fully charged, the lights on the device will turn off.

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