Product name: 7days Medicine Box
Item No: 101870 v001: 20/08/2025

Please read and retain these instructions for future reference

3PAGEN

GB INSTRUCTION MANUAL

Please read and retain these instructions for future reference

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INSTRUCTIONS FOR USE:

STEP 1: Prepare Your Medications

- Wash and dry your hands before handling pills.
- Check your prescription/dosage instructions.

STEP 2: Fill the Organizer

Gently pull out the drawer for the desired day (e.g., "Montag") from the bottom.

Load Your Pills.

Place the correct dosage into each of the 4 compartments:

AM (Morgen)

NOON (Mittag)

PM (Abend)

BEDTIME (Nacht)

Close the Drawer Securely

Push the drawer back to the first layer.

STEP 3: Taking Your Medication

- Each morning, pull out the current day's drawer.
- Open only the needed compartment (e.g., "AM" for morning pills).
- Close the compartment immediately after use to maintain moisture protection.

STEP 4: Weekly Reset & Cleaning

- After Sunday, refill all drawers for the next week.
- Clean with damp cloth & mild soap. Do not submerge in water.
- Air dry completely before refilling.

CARE AND CLEANING INSTRUCTIONS:

- Clean the organizer regularly to prevent residue buildup. Wash each compartment and lid with mild soap and lukewarm water. Use a soft brush if necessary.

DISPOSAL

- 3PAGEN use recyclable or recycled packaging where possible.
- Please dispose of all packaging, paper, cartons in accordance with your local authority recycling regulations.
- At the end of the product's lifespan please check with your local council authorised household waste recycling centre for disposal.

Code & Symbol			
Type of Plastic	Polypropylene	Polystyrene	ABS

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