Product name: Himalaya Salt Wrap

Item No: 729190

v001: 11/02/2026

Please read and retain these instructions for future reference

3PAGEN

GB Operating Instructions:

Neck and Shoulder Himalayan Salt Wrap HOT & COLD THERAPY

Heat Therapy:

1. Place the Himalayan Salt Wrap in the microwave.
2. Heat according to your microwave's wattage as follows:
 1. 500 to 600 watts: 2 minutes 30 seconds
 2. 600 to 800 watts: 2 minutes 15 seconds
 3. 800 to 900 watts: 2 minutes
 4. 900 to 1000 watts: 1 minute 45 seconds
3. After heating, remove the Himalayan Salt Wrap from the microwave and check the temperature before use.
4. Wear the Himalayan Salt Wrap.

Cold Therapy:

- For cold therapy, place the Himalayan Salt Wrap in a plastic bag and store it in the freezer for at least one hour before use.

Important Notes:

- Never leave the Himalayan Salt Wrap unattended in the microwave.
- Heating times may vary depending on the microwave's wattage and brand.
- Overheating the Himalayan Salt Wrap can cause a fire hazard.

Caution:

- Always check the wrap's temperature before applying it to your body to ensure safety.
- Keep and review these instructions for future use.
- Allow the Himalayan Salt Wrap to fully cool before reheating.
- Do not exceed the recommended maximum heating time.
- This product is not a toy. Adult supervision is advised. Not suitable for children under 3 years old or individuals with medical conditions.
- May cause allergic skin reactions. Wash your skin thoroughly with water after handling.
- Always wash your hands after use.
- If you experience an allergic reaction, seek immediate medical attention.
- Do not use the wrap for heating bedding, or duvets, or as a substitute for a hot water bottle.

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