Rechargeable EMS Foot Massager

Please read and retain these instructions for future reference

Art: 75006

Instructions for Use

Controller x 1pc
USB Cable x 1pc
Massage mat x 1pc

The massage mat is not to be used by anyone with the following health conditions:

Anyone fitted with a pace make or automatic implanted cardiac defibrillator(AICD)

Pregnant

Suffering with deep vein thrombosis(DVT)

Anyone suffering with high blood pressure or heart disease

Plun in the controller through USB charging cable (Included)
Connect the USB charging cable with AC/DC adaptor(Not included)
The LED indicator light on the controller will be on twinkling

Press the “ + ” button of the cortroller to intensity power levels and power on (15 Levels)

Press the “ - ” button of the cortroller to decrease power levels and power off

Press the “ M ” button of the controller to choose your modes (6 Modes)

Automatically turns off after 15 seconds (If you don't have your feet on the mat)

Rechargeable EMS Foot Massager

Please read and retain these instructions for future reference

Art: 75006

Instructions for Use

BATTERIES

- **ALWAYS USE ALKALINE BATTERIES. DO NOT USE LITHIUM BATTERIES.**
- INSTALL BATTERIES IN ACCORDANCE WITH ITS POLARITY SYMBOLS. DO NOT FIT BATTERIES IN REVERSE.
- ENSURE THE BATTERY COMPARTMENT IS RE-SECURED BY THE SCREW OR SECURING DEVICE.
- DO NOT SHORT CIRCUIT THE BATTERY TERMINALS.
- REMOVE THE BATTERIES IF LEFT FOR LONG PERIODS.
- **KEEP BATTERIES AWAY FROM CHILDREN AND PETS.**
- IF YOU THINK BATTERIES MIGHT HAVE BEEN SWALLOWED OR PLACED INSIDE ANY PART OF THE BODY, SEEK IMMEDIATE MEDICAL ATTENTION.
- DO NOT MIX USED AND NEW BATTERIES.
- **NEVER DISPOSE OF BATTERIES IN A FIRE.**
- BATTERY DISPOSAL, SPENT OR EXPIRED BATTERIES MUST BE DISPOSED OF AND RECYCLED IN COMPLIANCE WITH YOUR LOCAL AUTHORITY REGULATIONS.
- FOLLOW THE BATTERY MANUFACTURER'S SAFETY, USAGE, AND DISPOSAL INSTRUCTIONS.

CARE AND CLEANING INSTRUCTIONS

- ALWAYS SWITCH OFF AND REMOVE BATTERIES BEFORE CLEANING.
- DO NOT IMMERSE IN WATER.
- DO NOT USE ABRASIVE CLEANERS OR SHARP UTENSILS TO CLEAN THIS APPLIANCE.
- WIPE SURFACES WITH A SOFT DRY CLOTH. DRY THOROUGHLY BEFORE FURTHER USE.
- STORE THE APPLIANCE IN A COOL, DRY PLACE.

DISPOSAL

- Coopers of Stortford use recyclable or recycled packaging where possible.
- Please dispose of all packaging, paper, cartons in accordance with your local authority recycling regulations.
- At the end of the product's lifespan please check with your local council authorised household waste recycling centre for disposal.
- Plastics, polybags – this is made from the following recyclable plastic.

Code & Symbol	 LDPE	 OTHER
Type of Plastic	Low Density Polyethylene	ABS / PU / POLYESTER

UK
CA

CE

CUSTOMER SERVICES:

Customers in GB (England, Scotland and Wales), please contact our customer services - UK Head Office, 11 Bridge Street, Bishop's Stortford, CM23 2JU, Tel: 0330 331 0300
Customers in Northern Ireland and Europe, please contact our customer services - D.S.B. LD, 160 Bd. de Fourmies, 59100 Roubaix, France, Tel: (+44) 1279 701269

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